

E R R A T U M

In our publication

N. Niederstebruch & D. Sixt

Standard Nutrient Agar 1 as a substitute
for blood-supplemented Müller-Hinton agar
for antibiograms in developing countries

Eur J Clin Microbiol Infect Dis
DOI 10.1007/s10096-012-1735-2

Received: 24 July 2012 / Accepted: 14 August 2012

on page 2, paragraph 1

the sentence

The blood has to be blended with the
agar matrix at temperatures below **40** °C, as blood will
denaturate at higher temperatures.

must be

The blood has to be blended with the
agar matrix at temperatures below **50** °C, as blood will
denaturate at higher temperatures.

- **signed in the original letter** -

Dipl.-Ing. Nils Niederstebruch

Chairman GLOBOLAB e.V.